

THE WORLD'S MOST UNIQUE BEAUTY MAGAZINE

NEWBEAUTY

LightStim
Featured in This
Issue of
NewBeauty

THE NUMBER ONE
PRODUCT
THAT CONQUERS
AGING

page 62

SCULPT
THE PERFECT
BODY

MELT FAT IN MINUTES

page 148

* LOOK 10 YEARS
YOUNGER
IN 2010

THE NEWEST ANTI-AGING BREAKTHROUGHS page 178





newbeauty.com

THE TECHNOLOGY COLLAGEN-STIMULATING LED LIGHTS

HOW IT WORKS Unlike most at-home LED devices, LightStim's Anti-Aging Light uses noninvasive amber- and red-colored light and infrared LED lights or infrared wavelengths to penetrate deep within the dermal layers of skin to engage fibroblast cells. Once these cells are stimulated, collagen and elastin production is excited, pore size is reduced, fine lines and wrinkles are less visible and skin tone and texture improve. "When the lights stimulate the cells, they begin to produce collagen and elastin, the dermis starts to slowly plump up, and folds, lines and wrinkles become less visible," says Steve Marchese, president of LightStim. **\$349, lightstim.com**

USE IT At least five to six times a week until you achieve the desired result, and then two to three times a week for maintenance.


“**WHAT THE EXPERT SAYS** As long as the intensity and strength of the red light is powerful enough, it can help heal the skin, combat pigmentation and even create a bit of collagen at the surface,” says Somerville.”



After eight weeks of use, fine lines and wrinkles are less noticeable. Images provided by LightStim.

"This handheld covers a bigger surface area and is easier to handle than other similar devices we've used. We definitely noticed some minimizing of wrinkles, but some users might find it difficult to work it into their daily routine."

the handheld laser
LIGHTSTIM



the instant eye brightener
CLARISONIC OPAL SONIC INFUSION SYSTEM

"So easy to use! We like how the serum is actually in the tool, so it eliminates the extra step of applying a product. The area under our eyes was instantly brighter and it felt firmer too. Two hours after using the device, the skin around our eyes looked more 'awake' than usual. The only drawback is that the vibrations of the device are intense."

“**WHAT THE EXPERT SAYS** For patients seeking an at-home device that reduces wrinkles and fine lines, I recommend the Opal Sonic Infusion System,” says Dr. Shamban. “When you deliver nutritional ingredients deep into the skin, you see quicker results that last longer.”

THE TECHNOLOGY SONIC MICROMASSAGE

HOW IT WORKS From the innovators of sonic technology comes the new Opal Sonic Infusion System, which uses sonic micromassage to wake up tired-looking eyes. To use, apply the tool to the eye area—a marine- and botanical-rich serum is already infused into the device—and move it in a circular motion around the bone of the eye. According to the company, when the system is used under each eye for 20 seconds, the eyes appear brighter, refreshed and hydrated, making fine lines and wrinkles less apparent. **\$245, clarisonic.com**

USE IT Twice a day